

Improved pulmonary resistance in health sheep following Targeted Lung Denervation (TLD)



ML Mayse¹, J Hummel², PJ Johnson¹, K.T. Rouw¹
¹Nuvaira, Inc., Minneapolis, MN, USA
²University of North Carolina, Chapel Hill, NC, USA

Introduction: Targeted lung denervation (TLD) is a bronchoscopic procedure developed to relieve obstructive airways disease by disrupting pulmonary parasympathetic inputs to the lung. A denervation system and novel radiofrequency ablation catheter with dual cooled design was developed to facilitate ablation of bronchial branches of the vagal trunk while minimizing airway wall effects.

Aim: To determine the effect of lung denervation on pulmonary resistance in healthy sheep.

Methods: Following IACUC approval, 3 sheep underwent circumferential ablation of both main bronchi using a lung denervation system (Nuvaira™, Inc., USA). All sheep were mixed breed and weighed 47-51 Kg. Prior to and following TLD therapy, each animal underwent a series of pulmonary resistance measurements using a custom built forced oscillometry system. Atropine was utilized as a positive control to demonstrate bronchodilation prior to TLD.

Results: The effect of TLD therapy was evident in all three animals. On average, the airway resistance decreased by 27% (1.34 cmH₂O/(L/s)) and 30% (1.28 cmH₂O/(L/s)) in response to the peak effect of atropine and post-TLD therapy, respectively. The effect from TLD was maintained at 90 min post treatment while the effect of atropine resolved over the same time frame.


Conclusion: Healthy sheep subjected to TLD demonstrated decreased airway resistance comparable to the peak effect of atropine.

Background – Targeted lung denervation (TLD) utilizes radiofrequency (RF) ablation within the main bronchi to disrupt nerve inputs to the lungs

- Parasympathetic motor nerve tone in COPD
 - Patients with COPD have increased parasympathetic tone in the lungs¹
 - Anticholinergics decrease parasympathetic tone and reduce resistance to airflow in patients with COPD^{2,3}
- Targeted lung denervation (TLD)
 - A novel one time anticholinergic bronchoscopy procedure developed by Nuvaira™ Inc.
 - Utilizes circumferential RF ablation within the main bronchi to disrupt parasympathetic motor nerve input to the lungs

– RF catheter is designed to achieve nerve disruption with minimal effects to the airway surface and peribronchial structures

- The Nuvaira™ Dual Cooled RF Catheter:
 - Provides conductive cooling to the airway surface during RF energy delivery
 - Maintains the integrity of the airway mucosa during TLD
 - Targets bronchial nerve fascicles at depth from the airway surface
- TLD using the Nuvaira™ Catheter is designed to minimize ablation effects to non targeted surrounding lung tissue




References –
¹Undem B. and Kollarik M. The Role of Vagal Afferent Nerves in Chronic Obstructive Pulmonary Disease. Proc Am Thor Soc 2005; 2: 355-360.
²Belmonte K. Cholinergic Pathways in the Lungs and Anticholinergic Therapy for Chronic Obstructive Pulmonary Disease. Proc Am Thor Soc 2005; 2: 297-305.
³Qaseem A, Wilt T, Weinberger S, et al. Diagnosis and management of stable chronic obstructive pulmonary disease: A clinical practice Guideline update from the ACP, ACCP, ATS and ERS. Ann Intern Med. 2011; 155: 179-91.

Hypothesis – TLD with the Nuvaira™ Lung Denervation System in sheep disrupts parasympathetic input to the lungs decreasing pulmonary resistance.


Methods – 3 sheep underwent circumferential ablation of both main bronchi using a lung denervation system (Nuvaira™, Inc., USA).

- Each bronchus of 3 sheep underwent TLD with the Nuvaira™ System
 - The study was conducted under the guidance of an Institutional Animal Care and Use Committee in accordance with the study facility SOP and Animal Welfare Act of 1966
 - Each bronchus (n=6) underwent 8 sequential RF activations evenly spaced around the circumference of the airway (octants)
 - Prior to and following TLD therapy, each animal underwent a series of pulmonary resistance measurements using a custom built forced oscillometry system.
 - Atropine was utilized as a positive control to demonstrate bronchodilation prior to TLD.


Pulmonary Resistance Measurement




- A forced oscillation system to measure the pulmonary resistance was assembled
- Breathing of the animal was controlled by ventilator set to 10-15 bpm, VT 0.5 L, P_{insp} = 10 cm H₂O, and PEEP = 5cm H₂O
- A 5Hz pressure wave was generated by a custom speaker box and superimposed on the animals' controlled breathing.
- Air flow was measured by a calibrated, heated pneumotachograph, and pressure was measured by a pressure transducer
- Pneumotachograph and pressure transducer output were displayed on an oscilloscope, 5 Hz component isolated, and airway impedance determined by dividing airflow by pressure



Results – Atropine caused a transient 26% peak reduction in pulmonary resistance that normalized after 25min



TLD using the lung denervation system caused a sustained >30% reduction in pulmonary resistance in treated animals



Conclusion – TLD treatment induces bronchodilation through nerve disruption that is sustained and comparable in magnitude to the anticholinergic effect produced by atropine.